

real-time, D.I.Y., laboratory grade potency data to any sector within the cannabis ecosystem. Now you can get highly accurate, repeatable results in a flash, at a fraction of the cost.

TECH SPECS

PRICING


BROCHURE DOWNLOAD



The Beacon - winner of the Award



## BEACON™ TECHNICAL SPECIFICATIONS



Beacon	
<b>Cannabinoid Measurements</b>	Total potential THC, THC-A, Delta 9 THC Total potential CBD and CBD-A
<b>Dimensions (L x W x H)</b>	220 x 160 x 210 mm
<b>Weight (lbs)</b>	4.2
<b>Power Consumption (W)</b>	15
<b>Input Voltage (VDC)</b>	12
<b>Wavelength Ranges</b>	NIR (1550-1985nm)
<b>Lamp Type</b>	Tungsten Halogen
<b>Lamp Life (Hrs)</b>	2000
<b>User Interface</b>	On-Board Touchscreen
<b>Communications</b>	Ethernet, WiFi, USB
<b>Construction Materials</b>	Powdercoated Aluminum Window is sapphire



- Our Science
- Profiler
- Extract Producers
- Blog
- info@sageanalytics.com
- Our Service Promise
- Laboratories
- Articles
- Legal
- Dispensaries
- Videos

