

housing was designed to be impact resistant and even power washed for easy clean-up.

TECH SPECS

PRICING


BROCHURE

X



LUMINARY™ PROFILER

TECHNICAL SPECIFICATIONS



Luminary™ Profiler	
Cannabinoid Measurements	Total potential THC, THC-A, Delta 9 THC Total potential CBD and CBD-A
Dimensions (L x W x H)	230 x 340 x 160 mm
Weight (lbs)	11.5
Power Consumption (W)	15
Input Voltage (VDC)	12
Wavelength Ranges	NIR (1550–1985nm)
Lamp Type	Tungsten Halogen
Lamp Life (Hrs)	2000
User Interface	iPad Mini
Communications	Ethernet, WiFi, USB
Construction Materials	Top/Bottom Plates-Anodized Aluminum Center Section-Cast Urethane Window is Sapphire



Our Story

Beacon

Growers

News

720-282-4350

Our Science

Profiler

Extract Producers

Blog

info@sageanalytics.com

Our Service Promise

Laboratories

Articles

Legal

Dispensaries

Videos